

## Atlas 2 mM Nucleotide Mix - Data Sheet

Cat. No.	Pack Size	Conc.
BHS0002	100 µl SAMPLE	2 mM each
BH40202	1 ml	2 mM each

Lot no:

Exp. Date:

Storage: -20 °C

### Applications:

- Standard PCR
- qPCR
- cDNA synthesis
- DNA sequencing
- DNA labeling

### Description:

This **ready-to-use** Atlas Nucleotide Mix consists of dATP, dCTP, dGTP and dTTP (each at a final concentration of 2 mM) in sterile deionized water, within one tube.

### Quality Control

- Functionally tested in PCR with Bioatlas DNA Polymerases.
- Purity of each component is  $\geq 99\%$  (HPLC).
- Free of endo-, exodeoxyribonuclease, ribo-nuclease, phosphatase and nicking activities.

### Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only.

*Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.*